

For Research Use Only

ACY3 Monoclonal antibody

Catalog Number: 68711-1-Ig



Basic Information

Catalog Number: 68711-1-Ig	GenBank Accession Number: BC008689	Purification Method: Protein G purification
Size: 1000 µg/ml	GeneID (NCBI): 91703	CloneNo.: 1D7E2
Source: Mouse	UNIPROT ID: Q96HD9	Recommended Dilutions: WB 1:2000-1:10000
Isotype: IgG1	Full Name: aspartoacylase (aminocyclase) 3	
Immunogen Catalog Number: AG10234	Calculated MW: 319 aa, 35 kDa	
	Observed MW: 35 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : Calu-3 cells, A549 cells, Caco-2 cells
Species Specificity: Human	

Background Information

ACY3 is a 319 amino acid protein that deacetylates mercapturic acids in the proximal tubules of the kidney, where it is predominantly expressed.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68711-1-Ig (ACY3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.