

For Research Use Only

ACY3 Monoclonal antibody, PBS Only



Catalog Number: 68711-1-PBS

Basic Information

Catalog Number: 68711-1-PBS	GenBank Accession Number: BC008689	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 91703	CloneNo.: 1D7E2
Source: Mouse	UNIPROT ID: Q96HD9	
Isotype: IgG1	Full Name: aspartoacylase (aminocyclase) 3	
Immunogen Catalog Number: AG10234	Calculated MW: 319 aa, 35 kDa	
	Observed MW: 35 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Background Information

ACY3 is a 319 amino acid protein that deacetylates mercapturic acids in the proximal tubules of the kidney, where it is predominantly expressed.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68711-1-Ig (ACY3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68711-1-PBS in a different storage buffer formulation.