

For Research Use Only

# HMMR Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68718-2-PBS



## Basic Information

Catalog Number:

68718-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG33723

GenBank Accession Number:

NM\_012484

GeneID (NCBI):

3161

UNIPROT ID:

O75330

Full Name:

hyaluronan-mediated motility  
receptor (RHAMM)

Calculated MW:

84 kDa

Purification Method:

Protein A purification

CloneNo.:

1D12B3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50514-1, HMMR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68718-2-PBS. Detection antibody: 68718-3-PBS. Standard: Ag33723. Range: 0.391-25 ng/mL