

For Research Use Only

LDB1 Monoclonal antibody, PBS Only



Catalog Number: 68721-1-PBS

Basic Information

Catalog Number:

68721-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG34563

GenBank Accession Number:

BC000482

GeneID (NCBI):

8861

UNIPROT ID:

Q86U70

Full Name:

LIM domain binding 1

Calculated MW:

46 kDa

Purification Method:

Protein G purification

CloneNo.:

2F11C3

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68721-1-Ig (LDB1 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68721-1-PBS in a different storage buffer formulation.