

For Research Use Only

C15orf24 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68722-3-PBS



Basic Information

Catalog Number:

68722-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG26813

GenBank Accession Number:

BC104934

GeneID (NCBI):

56851

UNIPROT ID:

Q9NPA0

Full Name:

chromosome 15 open reading frame

24

Purification Method:

Protein G Magarose purification

CloneNo.:

1G7G8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50844-1, C15orf24 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68722-3-PBS. Detection antibody: 68722-4-PBS. Standard: Ag26813. Range: 0.391-100 ng/mL.