

For Research Use Only

ADAM17 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68725-2-PBS



Basic Information

Catalog Number: 68725-2-PBS	GenBank Accession Number: BC136783	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6868	CloneNo.: 1D5A9
Source: Mouse	UNIPROT ID: P78536	
Isotype: IgG2a	Full Name: ADAM metallopeptidase domain 17	
Immunogen Catalog Number: AG34022	Calculated MW: 824 aa, 93 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP50392-1, ADAM17 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68725-2-PBS. Detection antibody: 68725-3-PBS. Standard:Ag34022. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50392-2, ADAM17 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68725-2-PBS. Detection antibody: 68725-1-PBS. Standard:Ag34022. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50392-3, ADAM17 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68725-2-PBS. Detection antibody: 68725-4-PBS. Standard:Ag34022. Range: 0.098-100 ng/mL.