

For Research Use Only

# SYT15 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68729-2-PBS

## Basic Information

Catalog Number:

68729-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27264

GenBank Accession Number:

BC139914

GeneID (NCBI):

83849

UNIPROT ID:

Q9BQS2

Full Name:

synaptotagmin XV

Calculated MW:

421 aa, 47 kDa

Purification Method:

Protein G purification

CloneNo.:

1C8B5

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50014-1, SYT15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68729-1-PBS. Detection antibody: HRP-conjugated 68729-2-PBS. Standard: Ag27264. Range: 0.156-20 ng/mL.