

For Research Use Only

SAMHD1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68732-1-PBS

Basic Information

Catalog Number:

68732-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG18766

GenBank Accession Number:

BC036450

GeneID (NCBI):

25939

UNIPROT ID:

Q9Y3Z3

Full Name:

SAM domain and HD domain 1

Calculated MW:

626 aa, 72 kDa

Purification Method:

Protein G purification

CloneNo.:

1D10H7

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50020-1, SAMHD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68732-1-PBS. Detection antibody: HRP-conjugated 68732-2-PBS. Standard: Ag18766. Range: 0.313-20 ng/mL.