

For Research Use Only

# IFN beta Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68734-2-PBS

## Basic Information

Catalog Number:

68734-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

HZ-1298

GenBank Accession Number:

GeneID (NCBI):

3456

Full Name:

IFN beta

Purification Method:

Protein G purification

CloneNo.:

1B3F5

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50087-1, IFN beta Monoclonal Matched Antibody Pair, PBS Only.  
Capture antibody: 68734-1-PBS. Detection antibody: HRP-conjugated 68734-2-PBS.  
Standard:HZ-1298. Range: 1-50ng/mL