

For Research Use Only

Medin Monoclonal antibody

Catalog Number: 68736-1-Ig



Basic Information

Catalog Number: 68736-1-Ig	GenBank Accession Number: NM_005928	Purification Method: Protein A purification
Size: 1000 µg/ml	GeneID (NCBI): 4240	CloneNo.: 1H3D7
Source: Mouse	UNIPROT ID: Q08431	Recommended Dilutions: WB 1:2000-1:10000
Isotype: IgG1	Full Name: milk fat globule-EGF factor 8 protein	
Immunogen Catalog Number: AG30145	Calculated MW: 43 kDa	
	Observed MW: 51 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : human placenta tissue,
Species Specificity: Human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68736-1-Ig (MFG8 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.