

For Research Use Only

IL5 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68738-1-PBS

Basic Information

Catalog Number:

68738-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG11840

GenBank Accession Number:

BC066282

GeneID (NCBI):

3567

UNIPROT ID:

P05113

Full Name:

interleukin 5 (colony-stimulating factor, eosinophil)

Calculated MW:

134 aa, 15 kDa

Purification Method:

Protein G purification

CloneNo.:

1B8A8

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50091-1, IL5
Monoclonal Matched Antibody Pair, PBS Only.
Capture antibody: 68737-1-PBS. Detection
antibody: HRP-conjugated 68738-1-PBS.
Standard:HZ-1324. Range: 0.977-250 ng/mL