

For Research Use Only

# ALOX15B Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68743-1-PBS

## Basic Information

Catalog Number: 68743-1-PBS	GenBank Accession Number: BC035217	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 247	CloneNo.: 1A8C4
Source: Mouse	UNIPROT ID: O15296	
Isotype: IgG2a	Full Name: arachidonate 15-lipoxygenase, type B	
Immunogen Catalog Number: AG30796	Calculated MW: 676 aa, 76 kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP50096-1, ALOX15B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68743-1-PBS. Detection antibody: 68743-2-PBS. Standard: Ag30796. Range: 0.049-6.25 ng/mL