

For Research Use Only

ALOX15B Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68743-2-PBS

Basic Information

Catalog Number:

68743-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG30796

GenBank Accession Number:

BC035217

GeneID (NCBI):

247

UNIPROT ID:

O15296

Full Name:

arachidonate 15-lipoxygenase, type B

Calculated MW:

676 aa, 76 kDa

Purification Method:

Protein G purification

CloneNo.:

1A2B9

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50096-1, ALOX15B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68743-1-PBS. Detection antibody: 68743-2-PBS. Standard: Ag30796. Range: 0.049-6.25 ng/mL