

For Research Use Only

# PD-L1/CD274 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68745-2-PBS

## Basic Information

Catalog Number:

68745-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC074984

GeneID (NCBI):

29126

UNIPROT ID:

Q9NZQ7

Full Name:

CD274 molecule

Calculated MW:

290 aa, 33 kDa

Purification Method:

Protein G purification

CloneNo.:

2D4A12

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50090-1, PD-L1/CD274 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68745-1-PBS. Detection antibody: HRP-conjugated 68745-2-PBS. Standard:Eg31414. Range:0.049-1.563 ng/mL.