

For Research Use Only

CD80/B7-1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68746-1-PBS

Basic Information

Catalog Number: 68746-1-PBS	GenBank Accession Number: NM_005191	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 941	CloneNo.: 1A10C6
Source: Mouse	UNIPROT ID: P33681	
Isotype: IgG1	Full Name: CD80 molecule	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50098-1, CD80/B7-1 Monoclonal Matched Antibody Pair - BSA and Azide free. Capture antibody: 68746-1-PBS. Detection antibody: 68746-2-PBS. Standard:Eg0060. Range: 0.391-100 ng/mL