

For Research Use Only

TRIM38 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68747-1-PBS

Basic Information

Catalog Number: 68747-1-PBS	GenBank Accession Number: BC026930	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 10475	CloneNo.: 2B2D9
Source: Mouse	UNIPROT ID: O00635	
Isotype: IgG2b	Full Name: tripartite motif-containing 38	
Immunogen Catalog Number: AG25204	Calculated MW: 465 aa, 53 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50100-1, TRIM38 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68747-1-PBS. Detection antibody: 68747-2-PBS. Standard: Ag25204. Range: 1.563-25 ng/mL