

For Research Use Only

CTLA4 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68751-2-PBS

Basic Information

Catalog Number:

68751-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

BC070162

GeneID (NCBI):

1493

UNIPROT ID:

P16410

Full Name:

cytotoxic T-lymphocyte-associated protein 4

Calculated MW:

223 aa, 25 kDa

Purification Method:

Protein A purification

CloneNo.:

1G2C12

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50106-2, CTLA4 Monoclonal Matched Antibody Pair - BSA and Azide free. Capture antibody: 68751-1-PBS. Detection antibody: 68751-2-PBS. Standard:Eg0018. Range: 0.391-50 ng/mL.