

For Research Use Only

# TIGIT Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68752-2-PBS

## Basic Information

<b>Catalog Number:</b> 68752-2-PBS	<b>GenBank Accession Number:</b> NM_173799.4	<b>Purification Method:</b> Protein G purification
<b>Size:</b> 1mg/ml	<b>GeneID (NCBI):</b> 201633	<b>CloneNo.:</b> 2D10G6
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q495A1-1	
<b>Isotype:</b> IgG1	<b>Full Name:</b> T cell immunoreceptor with Ig and ITIM domains	
	<b>Calculated MW:</b> 26 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50107-1, TIGIT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68752-1-PBS. Detection antibody: 68752-2-PBS. Standard:Eg0748. Range: 3.125-100 ng/mL.