

For Research Use Only

OAT Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68753-1-PBS

Basic Information

Catalog Number:

68753-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG34783

GenBank Accession Number:

BC000964

GeneID (NCBI):

4942

UNIPROT ID:

P04181

Full Name:

ornithine aminotransferase (gyrate atrophy)

Calculated MW:

62 kDa

Purification Method:

Protein G purification

CloneNo.:

1G2D5

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50109-1, OAT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68753-1-PBS. Detection antibody: 68753-2-PBS. Standard: Ag34783. Range: 0.098-100 ng/mL.