

For Research Use Only

# LEP Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68756-2-PBS

## Basic Information

Catalog Number:

68756-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC060830

GeneID (NCBI):

3952

UNIPROT ID:

P41159

Full Name:

leptin

Calculated MW:

167 aa, 19 kDa

Purification Method:

Protein G purification

CloneNo.:

1E9A6

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50112-1, LEP Monoclonal Matched Antibody Pair - BSA and Azide free. Capture antibody: 68756-1-PBS. Detection antibody: 68756-2-PBS. Standard: Eg0328. Range: 0.098-100 ng/mL