

For Research Use Only

Leptin Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68756-5-PBS



Basic Information

Catalog Number:

68756-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC060830

GeneID (NCBI):

3952

UNIPROT ID:

P41159

Full Name:

leptin

Calculated MW:

167 aa, 19 kDa

Purification Method:

Protein G purification

CloneNo.:

1G3G6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50112-3, Leptin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68756-1-PBS. Detection antibody: 68756-5-PBS. Standard:Eg0328. Range: 0.195-100 ng/mL.



Cytometric bead array standard curve of MP50112-4, Leptin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68756-6-PBS. Detection antibody: 68756-5-PBS. Standard:Eg0328. Range: 0.781-50 ng/mL.