

For Research Use Only

IL1A Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68758-1-PBS

Basic Information

Catalog Number:

68758-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

HZ-1320

GenBank Accession Number:

BC013142

GeneID (NCBI):

3552

UNIPROT ID:

P01583

Full Name:

interleukin 1, alpha

Calculated MW:

271 aa, 31 kDa

Purification Method:

Protein G purification

CloneNo.:

1B9F5

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50113-1, IL1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68757-1-PBS. Detection antibody: 68758-1-PBS. Standard:HZ-1320. Range: 0.098-3.125 ng/mL.