

For Research Use Only

EPHA2 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68759-2-PBS

Basic Information

Catalog Number:

68759-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC037166

GeneID (NCBI):

1969

UNIPROT ID:

P29317

Full Name:

EPH receptor A2

Calculated MW:

976 aa, 108 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2C5E5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50114-1, EPHA2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68759-1-PBS. Detection antibody: 68759-2-PBS. Standard:Eg0146. Range: 0.195-50 ng/mL