

For Research Use Only

VEGFR1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68762-1-PBS

Basic Information

Catalog Number:
68762-1-PBS

Size:
1mg/ml

Source:
Mouse

Isotype:
IgG1

GenBank Accession Number:
NM_002019

GeneID (NCBI):
2321

UNIPROT ID:
P17948-1

Full Name:
fms-related tyrosine kinase 1
(vascular endothelial growth
factor/vascular permeability factor
receptor)

Calculated MW:
151kd

Purification Method:
Protein G purification

CloneNo.:
1G7C3

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50118-1, VEGFR1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68762-1-PBS. Detection antibody: 68762-2-PBS. Standard: Eg0438. Range: 6.25-100 ng/mL