

For Research Use Only

ANO1,DOG1 Monoclonal antibody, PBS Only (Detector)



Catalog Number:68764-2-PBS

Basic Information

Catalog Number:

68764-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG3320

GenBank Accession Number:

BC027590

GeneID (NCBI):

55107

UNIPROT ID:

Q5XXA6

Full Name:

anoctamin 1, calcium activated
chloride channel

Calculated MW:

986 aa, 114 kDa

Purification Method:

Protein A purification

CloneNo.:

1D3B2

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50123-1, ANO1, DOG1 Monoclonal Matched Antibody Pair - BSA and Azide free. Capture antibody: 68764-1-PBS. Detection antibody: 68764-2-PBS. Standard: Ag3320. Range: 0.098-100 ng/mL