For Research Use Only

EIF6 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68765-1-PBS

Basic Information

Catalog Number:

68765-1-PBS BC001119

 Size:
 GeneID (NCBI):

 1mg/ml
 3692

 Source:
 UNIPROT ID:

 Mouse
 P56537

Isotype: Full Name:

lgG1 eukaryotic translation initiation

GenBank Accession Number:

Immunogen Catalog Number: factor 6

AG30249 Calculated MW: 27 kDa

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information


Storage

Storage:

Store at -80°C. Storage Buffer: PBS Only Purification Method: Protein G purification

CloneNo.: 1A4B4

Selected Validation Data



Cytometric bead array standard curve of MP50124-1, EIF6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68765-1-PBS. Detection antibody: 68765-2-PBS. Standard:Ag30249. Range: 3.125-100 ng/mL