

For Research Use Only

EIF6 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68765-2-PBS

Basic Information

Catalog Number: 68765-2-PBS	GenBank Accession Number: BC001119	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 3692	CloneNo.: 1E2G5
Source: Mouse	UNIPROT ID: P56537	
Isotype: IgG1	Full Name: eukaryotic translation initiation factor 6	
Immunogen Catalog Number: AG30249	Calculated MW: 27 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50124-1, EIF6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68765-1-PBS. Detection antibody: 68765-2-PBS. Standard: Ag30249. Range: 3.125-100 ng/mL.