

For Research Use Only

EIF6 Monoclonal antibody, PBS Only

Catalog Number: 68765-3-PBS



Basic Information

Catalog Number:

68765-3-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG30249

GenBank Accession Number:

BC001119

GeneID (NCBI):

3692

UNIPROT ID:

P56537

Full Name:

eukaryotic translation initiation

factor 6

Calculated MW:

27 kDa

Purification Method:

Protein A purification

CloneNo.:

1F2A8

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68765-3-Ig (EIF6 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68765-3-PBS in a different storage buffer formulation.