

For Research Use Only

ANKRD13B Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68766-1-PBS

Basic Information

Catalog Number: 68766-1-PBS	GenBank Accession Number: BC032554	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 124930	CloneNo.: 1D2F11
Source: Mouse	UNIPROT ID: Q86YJ7	
Isotype: IgG2a	Full Name: ankyrin repeat domain 13B	
Immunogen Catalog Number: AG22725	Calculated MW: 626 aa, 70 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50125-1, ANKRD13B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68766-1-PBS. Detection antibody: 68766-2-PBS. Standard: Ag22725. Range: 1.563-100 ng/mL.