

For Research Use Only

GDF15 Monoclonal antibody, PBS Only (Capture/Detector)



Catalog Number: 68768-2-PBS

Basic Information

Catalog Number:

68768-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33750

GenBank Accession Number:

BC008962

GeneID (NCBI):

9518

UNIPROT ID:

Q99988

Full Name:

growth differentiation factor 15

Calculated MW:

34 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B5D12

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Sandwich ELISA standard curve of MP50133-10, GDF15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68768-6-PBS. Detection antibody: HRP-conjugated 68768-2-PBS. Standard:ag33750. Range: 1.563-100 ng/mL.




Sandwich ELISA standard curve of MP50133-11, GDF15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68768-2-PBS. Detection antibody: HRP-conjugated 68768-3-PBS. Standard:ag33750. Range: 1.563-100 ng/mL.



Sandwich ELISA standard curve of MP50133-5, GDF15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68768-1-PBS. Detection antibody: HRP-conjugated 68768-2-PBS. Standard:ag33750. Range: 1.563-100 ng/mL.



Sandwich ELISA standard curve of MP50133-6, GDF15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68768-3-PBS. Detection antibody: HRP-conjugated 68768-2-PBS. Standard:ag33750. Range: 1.563-100 ng/mL.



Sandwich ELISA standard curve of MP50133-7, GDF15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68768-4-PBS. Detection antibody: HRP-conjugated 68768-2-PBS. Standard:ag33750. Range: 1.563-100 ng/mL.