

For Research Use Only

TNFR2 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68771-2-PBS

Basic Information

Catalog Number: 68771-2-PBS	GenBank Accession Number: BC052977	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 7133	CloneNo.: 3D6E8
Source: Mouse	UNIPROT ID: P20333	
Isotype: IgG3	Full Name: tumor necrosis factor receptor superfamily, member 1B	
Immunogen Catalog Number: AG30497	Calculated MW: 48 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50139-1, TNFR2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68771-1-PBS. Detection antibody: HRP-conjugated 68771-2-PBS. Standard: Ag30497. Range: 0.098-25 ng/mL.