

For Research Use Only

# BMP-7 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68772-1-PBS

## Basic Information

Catalog Number:

68772-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

HZ-1229

GenBank Accession Number:

GeneID (NCBI):

655

Full Name:

bone morphogenetic protein 7

Purification Method:

Protein G purification

CloneNo.:

1A7C1

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50138-1, BMP-7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68772-1-PBS. Detection antibody: 68772-2-PBS. Standard:HZ-1229. Range: 3.125-100 ng/mL.