

For Research Use Only

BCMA/TNFRSF17 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68775-1-PBS



Basic Information

Catalog Number:

68775-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG26996

GenBank Accession Number:

BC058291

GeneID (NCBI):

608

UNIPROT ID:

Q02223

Full Name:

tumor necrosis factor receptor
superfamily, member 17

Calculated MW:

184 aa, 20 kDa

Purification Method:

Protein G purification

CloneNo.:

3F12C11

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50140-1, TNFRSF17 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68775-1-PBS. Detection antibody: HRP-conjugate 68775-2-PBS. Standard: Ag26996. Range: 0.049-25 ng/mL.