

For Research Use Only

# TNFRSF17 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68775-2-PBS

## Basic Information

<b>Catalog Number:</b> 68775-2-PBS	<b>GenBank Accession Number:</b> BC058291	<b>Purification Method:</b> Protein G purification
<b>Size:</b> 1mg/ml	<b>GeneID (NCBI):</b> 608	<b>CloneNo.:</b> 1B8C4
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q02223	
<b>Isotype:</b> IgG1	<b>Full Name:</b> tumor necrosis factor receptor superfamily, member 17	
<b>Immunogen Catalog Number:</b> AG26996	<b>Calculated MW:</b> 184 aa, 20 kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP50140-1, TNFRSF17 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68775-1-PBS. Detection antibody: HRP-conjugate 68775-2-PBS. Standard: Ag26996. Range: 0.049-25 ng/mL.