

For Research Use Only

PPP2R1A Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68778-1-PBS

Basic Information

Catalog Number:

68778-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8661

GenBank Accession Number:

BC001537

GeneID (NCBI):

5518

UNIPROT ID:

P30153

Full Name:

protein phosphatase 2 (formerly 2A),
regulatory subunit A, alpha isoform

Calculated MW:

589 aa, 65 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D1C1

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50144-1, PPP2R1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68778-1-PBS. Detection antibody: 68778-2-PBS. Standard: Ag8661. Range: 6.25-100 ng/mL