

For Research Use Only

HMGC51 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68779-5-PBS

Basic Information

| | | |
|---|---|--|
| Catalog Number: 68779-5-PBS | GenBank Accession Number: BC083514 | Purification Method: Protein G Magarose purification |
| Size: 1mg/ml | GeneID (NCBI): 3157 | CloneNo.: 2E5A1 |
| Source: Mouse | UNIPROT ID: Q01581 | |
| Isotype: IgG1 | Full Name: 3-hydroxy-3-methylglutaryl-Coenzyme A synthase 1 (soluble) | |
| Immunogen Catalog Number: AG11937 | Calculated MW: 520 aa, 57 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50145-3, HMGC51 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68779-4-PBS. Detection antibody: 68779-5-PBS. Standard: Ag11937. Range: 0.098-100 ng/mL.