

For Research Use Only

HMGC51 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68779-6-PBS

Basic Information

Catalog Number: 68779-6-PBS	GenBank Accession Number: BC083514	Purification Method: Protein G Magarose purification
Size: 1mg/ml	GeneID (NCBI): 3157	CloneNo.: 1E10A8
Source: Mouse	UNIPROT ID: Q01581	
Isotype: IgG1	Full Name: 3-hydroxy-3-methylglutaryl-Coenzyme A synthase 1 (soluble)	
Immunogen Catalog Number: AG11937	Calculated MW: 520 aa, 57 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50145-4, HMGC51 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68779-6-PBS. Detection antibody: 68779-2-PBS. Standard: Ag11937. Range: 0.098-100 ng/mL.