

For Research Use Only

PRPH2 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68780-2-PBS

Basic Information

Catalog Number: 68780-2-PBS	GenBank Accession Number: BC074720	Purification Method: Protein A Magarose purification
Size: 1mg/ml	GeneID (NCBI): 5961	CloneNo.: 1F6B3
Source: Mouse	UNIPROT ID: P23942	
Isotype: IgG2a	Full Name: peripherin 2 (retinal degeneration, slow)	
Immunogen Catalog Number: AG12555	Calculated MW: 346 aa, 39 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50146-1, PRPH2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68780-1-PBS. Detection antibody: 68780-2-PBS. Standard:Ag12555. Range: 1.563-100 ng/mL



Cytometric bead array standard curve of MP50146-2, PRPH2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68780-3-PBS. Detection antibody: 68780-2-PBS. Standard:Ag12555. Range: 0.781-100 ng/mL