

For Research Use Only

PRPH2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68780-3-PBS

Basic Information

| | | |
|---|--|--|
| Catalog Number: 68780-3-PBS | GenBank Accession Number: BC074720 | Purification Method: Protein A Magarose purification |
| Size: 1mg/ml | GeneID (NCBI): 5961 | CloneNo.: 3D6G10 |
| Source: Mouse | UNIPROT ID: P23942 | |
| Isotype: IgG2a | Full Name: peripherin 2 (retinal degeneration, slow) | |
| Immunogen Catalog Number: AG12555 | Calculated MW: 346 aa, 39 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50146-2, PRPH2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68780-3-PBS. Detection antibody: 68780-2-PBS. Standard: Ag12555. Range: 0.781-100 ng/mL.