

For Research Use Only

OX40L/TNFSF4 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68781-1-PBS

Basic Information

Catalog Number:

68781-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM_003326.5

GeneID (NCBI):

7292

UNIPROT ID:

P23510

Full Name:

tumor necrosis factor (ligand)
superfamily, member 4

Purification Method:

Protein G purification

CloneNo.:

1D4F2

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50147-1, OX40L/TNFSF4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68781-1-PBS. Detection antibody: 68781-2-PBS. Standard: Eg0608. Range: 1.563-100 ng/mL