

For Research Use Only

# OX40L/TNFSF4 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68781-2-PBS

## Basic Information

Catalog Number:

68781-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2b

GenBank Accession Number:

NM\_003326.5

GeneID (NCBI):

7292

UNIPROT ID:

P23510

Full Name:

tumor necrosis factor (ligand)  
superfamily, member 4

Purification Method:

Protein A purification

CloneNo.:

2C6A7

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50147-1, OX40L/TNFSF4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68781-1-PBS. Detection antibody: 68781-2-PBS. Standard: Eg0608. Range: 1.563-100 ng/mL