

For Research Use Only

PGDS Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68783-1-PBS

Basic Information

Catalog Number:

68783-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG34483

GenBank Accession Number:

BC005939

GeneID (NCBI):

5730

UNIPROT ID:

P41222

Full Name:

prostaglandin D2 synthase 21kDa
(brain)

Calculated MW:

21 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1F2B10

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50148-1, PGDS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68783-1-PBS. Detection antibody: 68783-2-PBS. Standard: Ag34483. Range: 3.125-100 ng/mL.