

For Research Use Only

PGDS Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68783-2-PBS

Basic Information

Catalog Number: 68783-2-PBS	GenBank Accession Number: BC005939	Purification Method: Protein A Magarose purification
Size: 1mg/ml	GeneID (NCBI): 5730	CloneNo.: 1C12E5
Source: Mouse	UNIPROT ID: P41222	
Isotype: IgG2a	Full Name: prostaglandin D2 synthase 21kDa (brain)	
Immunogen Catalog Number: AG34483	Calculated MW: 21 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50148-1, PGDS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68783-1-PBS. Detection antibody: 68783-2-PBS. Standard: Ag34483. Range: 3.125-100 ng/mL.