

For Research Use Only

ACTR8 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68784-2-PBS

Basic Information

Catalog Number:

68784-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG34088

GenBank Accession Number:

BC032744

GeneID (NCBI):

93973

UNIPROT ID:

Q9H981

Full Name:

ARP8 actin-related protein 8 homolog (yeast)

Calculated MW:

329 aa, 37 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3C4B11

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50150-1, ACTR8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68784-1-PBS. Detection antibody: 68784-2-PBS. Standard: Ag34088. Range: 0.781-50 ng/mL