

For Research Use Only

# RIPK3 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68786-3-PBS

## Basic Information

Catalog Number:

68786-3-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG11759

GenBank Accession Number:

BC062584

GeneID (NCBI):

11035

UNIPROT ID:

Q9Y572

Full Name:

receptor-interacting serine-threonine  
kinase 3

Calculated MW:

518 aa, 57 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3B6E7

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50152-2, RIPK3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68786-3-PBS. Detection antibody: 68786-4-PBS. Standard:ag11759. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP50152-3, RIPK3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68786-3-PBS. Detection antibody: 68786-2-PBS. Standard:ag11759. Range: 0.391-50 ng/mL.