

For Research Use Only

RNASEH1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68787-1-PBS

Basic Information

Catalog Number: 68787-1-PBS	GenBank Accession Number: BC002973	Purification Method: Protein G Magarose purification
Size: 1mg/ml	GeneID (NCBI): 246243	CloneNo.: 3C7F6
Source: Mouse	UNIPROT ID: O60930	
Isotype: IgG1	Full Name: ribonuclease H1	
Immunogen Catalog Number: AG33948	Calculated MW: 32 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50153-1, RNASEH1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68787-1-PBS. Detection antibody: 68787-2-PBS. Standard: Ag33948. Range: 0.391-100 ng/mL.