

For Research Use Only

# SLC35A1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68788-3-PBS

## Basic Information

Catalog Number:

68788-3-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG9524

GenBank Accession Number:

BC017807

GeneID (NCBI):

10559

UNIPROT ID:

P78382

Full Name:

solute carrier family 35 (CMP-sialic acid transporter), member A1

Calculated MW:

337 aa, 37 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2G12C2

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50154-2, SLC35A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68788-3-PBS. Detection antibody: 68788-2-PBS. Standard: Ag9524. Range: 1.563-50 ng/mL