

For Research Use Only

CEBPA Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68793-1-PBS

Basic Information

Catalog Number: 68793-1-PBS	GenBank Accession Number: NM_004364	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 1050	CloneNo.: 3C10G4
Source: Mouse	UNIPROT ID: P49715	
Isotype: IgG1	Full Name: CCAAT/enhancer binding protein (C/EBP), alpha	
Immunogen Catalog Number: AG34094	Calculated MW: 38 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50160-1, CEBPA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68793-1-PBS. Detection antibody: 68793-2-PBS. Standard: Ag34094. Range: 1.563-50 ng/mL