

For Research Use Only

# PD-L2/B7-DC Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68795-2-PBS

## Basic Information

Catalog Number:

68795-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC074766

GeneID (NCBI):

80380

UNIPROT ID:

Q9BQ51

Full Name:

programmed cell death 1 ligand 2

Calculated MW:

273 aa, 31 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3E2D7

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50162-1, PD-L2/B7-DC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68795-1-PBS. Detection antibody: 68795-2-PBS. Standard: Eg0507. Range: 0.781-100 ng/mL