

For Research Use Only

CD200R1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68796-1-PBS

Basic Information

Catalog Number:

68796-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC069743

GeneID (NCBI):

131450

UNIPROT ID:

Q8TD46

Full Name:

CD200 receptor 1

Purification Method:

Protein G Magarose purification

CloneNo.:

2C1B9

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50163-1, CD200R1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68796-1-PBS. Detection antibody: 68796-2-PBS. Standard:Eg0011. Range: 1.563-100 ng/mL.